



NICHIPET AUTOCLAVABLE

Thumb Knob Set (Cover, Knob & Screw)

G2-01G/01W/01Y/01B




Lock Handle (Lever)
5G2-0071




2nd Spring
5G2-520200
5G2-525000

2nd Spring Holder
5G2-530200
5G2-535000



Tool to remove the 2nd Spring Holder
5G2-JIG




Calibration Screw Assy
5G2-550000

Lock Screw
5G2-008





Volume Counter Assy
w/Stroke Pipe

G2-560200
G2-561000
G2-565000




Shaft
5G2-009


NICHIPET AUTOCLAVABLE



NP-5000/10ML



NP-2/10/20/100/200/1000



NICHIPET AUTOCLAVABLE